
# ProSeal™ In Line Pump Set

REF 423110

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
# ProSeal™ In Line Pump Set

REF 423110

Priming Volume 0.02 ml (Cannula Only)  
Gauge of Cannula 18G  
External Fitting Twist Off Port

Packaging  
50 pcs/Inner Box  
200 pcs/Outer Carton

NOT MADE with NATURAL RUBBER LATEX or DEHP

EC REP

MT Promedt Consulting GmbH  
Altenhofstr.80  
66386 St. Ingbert Germany  
+49-6894-581020

Epic International (Thailand) Co., Ltd.  
Hemaraj Eastern Seaboard Industrial Estate  
500/73 Moo 2, Tasit, Pluakdaeng, Rayong  
Thailand 21140

www.epic-med.com


## Indications for Use:

The ProSeal Closed System drug Transfer Device (CSTD) mechanically prohibits environmental contaminants from entering the system and the escape of drug or vapor concentrations from the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols, and spills. The ProSeal system also prevents the introduction of microbial contaminants into the drug or fluid path for up to 168 hours or 7 days.

Page 5 of 6

780062 Rev.01, Date of issue 2020-06-27

DIRECTIONS - Use Aseptic Technique

1. Attach \* In Line Pump Set to any Proseal Bag Spike.
2. Twist off the end cap of the port. Insert the spike of any IV administration set into the port. The In Line Pump Set could be detached and re-attached to any subsequent fluid or medication bag. In each step, a closed system fluid transfer is maintained.
3. Discard used In Line Pump Set without disconnecting from IV administration set in accordance to disposal procedures for biohazardous materials of your facility.

\* Ensure that the opening of In Line Pump Set is engaged with the end of spike before pushing to connect.

## Contraindications

The device is contraindicated whenever:

- The drug to be prepared is contraindicated to polycarbonate, ABS, stainless steel, polypropylene, PVC, silicone and polyisoprene.

## Warnings

- Use accepted IV and pharmacy practice.
- The performance of the self-sealing membrane of the device is reduced after multiple perforation s.
- Do not re-use to avoid contamination.

## Precautions

- Do not use when caps and /or components are loose.
- The device is sterile unless packaging is damaged. Do not use when packaging is damaged.

Page 6 of 6