
ProSeal™ In Line Pump Set

REF 423850

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
ProSeal™ In Line Pump Set

REF 423850

Priming Volume 0.02 ml (Cannula Only)
Gauge of Cannula 18G
External Fitting Spike Port
Compatible with DMA (Dimethylacetamide)

NOT MADE with NATURAL RUBBER LATEX or DEHP

MT Promedt Consulting GmbH
Altenhofstr.80
66386 St. Ingbert Germany
+49-6894-581020

Epic International (Thailand) Co., Ltd.
Hemaraj Eastern Seaboard Industrial Estate
500/73 Moo 2, Tasit, Pluakdaeng, Rayong
Thailand 21140

www.epic-med.com

Indications for Use:

The ProSeal Closed System drug Transfer Device (CSTD) mechanically prohibits environmental contaminants from entering the system and the escape of drug or vapor concentrations from the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols, and spills. The ProSeal system also prevents the introduction of microbial contaminants into the drug or fluid path for up to 168 hours or 7 days.

CE 0123

STERILE EO


Rx

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Page 5 of 8

780139 Rev.04, Date of issue 2024-01-06

DIRECTIONS - Use Aseptic Technique

1. Attach * In Line Pump Set to any Proseal Bag Spike.
2. Remove the spike port cover. Insert the spike of any IV administration set into the port. The In Line Pump Set could be detached and re-attached to any subsequent fluid or medication bag. In each step, a closed system fluid transfer is maintained.
3. Discard used In Line Pump Set without disconnecting from IV administration set in accordance to disposal procedures for biohazardous materials of your facility.

* Ensure that the opening of In Line Pump Set is engaged with the end of ProSeal Bag Spike before pushing to connect.

Contraindications

The device is contraindicated whenever:

- The drug to be prepared is contraindicated to copolyester, stainless steel, polypropylene, TPE and polyisoprene.

Warnings

- Use accepted IV and pharmacy practice.
- The performance of the self-sealing membrane of the device is reduced after multiple perforations.
- Do not re-use to avoid contamination.

Precautions

- Do not use when caps and /or components are loose.
- The device is sterile unless packaging is damaged. Do not use when packaging is damaged.

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