
## ToxiSeal Vial Adaptors with External Balloon and Valve

		1 pcs/pack	5 pcs/pack
Vial 13 mm	-	420300	420305
	Concave	420400	420405
Vial 20 mm	-	420310	420315
	Concave	420410	420415
Vial 28 mm	-	420320	420325
	Concave	420450	420455

External Balloon Pressure Equalization Volume 80 ml  
Air Filter Pore Size 0.1 µm  
Carbon Filter Activated Carbon Fabric  
Compatible with DMA (Dimethylacetamide)

NOT MADE with NATURAL RUBBER LATEX or DEHP

EC REP

MT Promed Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert Germany  
+49-6894-581020

Epic International (Thailand) Co., Ltd.  
Hemraj Eastern Seaboard Industrial Estate  
500/73 Moo 2, Tasit, Pluakdaeng, Rayong  
Thailand 21140

www.epic-med.com

### Indications for Use:

The ProSeal Closed System drug Transfer Device (CSTD) mechanically prohibits environmental contaminants from entering the system and the escape of drug or vapor concentrations from the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols, and spills. The ProSeal system also prevents the introduction of microbial contaminants into the drug or fluid path for up to 168 hours or 7 days.

CE 0123

STERILE EO

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Page 5 of 8

780121 Rev.03 - Date of issue: 2023-04-08

ToxiSeal Vial Adaptors seal the closure of the vial to which it is attached. The Air Vent assembly consists of 2 layers of filtration membranes, a hydrophobic 0.1 micron sterile filter and an activated carbon filter that absorbs drug aerosols and vapor.

The external balloon collects drug aerosols and volatile organic compounds released due to un-equalized pressure in the vial during compounding. The fluid path of ToxiSeal Vial Adaptors is made of Tritan™ Copolyester material, which is compatible to DMA (Dimethylacetamide) solvent in certain chemotherapy drugs.

The injection site on the ToxiSeal Vial Adaptor and all its corresponding interface membranes exhibit a dry connection with the communicating surfaces in a fluid transfer. The use of this component and its appropriate ProSeal CSTD connecting component reduces the risk of microbial ingress for up to 168 hours or 7 days.

### DIRECTIONS - Use Aseptic Technique

1. Remove the protective cap on vial spike and attach ToxiSeal Vial Adaptor firmly to vial.
2. Prior to every access, swab top of Injection Site with 70% isopropyl alcohol (15 seconds) and allow to dry (approximately 30 secs).  
Dry Time is dependent on temperature, humidity, ventilation area.
3. Attach ProSeal Injector or Injector Plus (Syringe Adaptor) to the Injection Site of ToxiSeal Vial Adaptor.
4. For drugs require reconstitution, hold the vial upright and inject diluent. Do not inject diluent more than the volume indicated for pressure equalization.
5. For liquid drugs, the one way valve allows solution to be withdrawn without injecting air into external balloon. Invert vial to withdraw fluid. Avoid pushing syringe plunger when the vial is inverted as excessive pressure may compromise the filter membranes.
6. Discard used ToxiSeal Vial Adaptor with vial intact in accordance to disposal procedures for biohazardous materials of your facility.

### Contraindications

The device is contraindicated whenever:

- The drug to be prepared is contraindicated to copolyester, stainless steel, polypropylene, TPE and polyisoprene.

### Warnings

- Use accepted IV and pharmacy practice.
- The performance of the self-sealing membrane of the device is reduced after multiple perforations.
- Do not re-use to avoid contamination.
- Only use the device model that matches the vial neck size and vial stopper type.

### Precautions

- Do not use when caps and/or components are loose.
- The device is sterile unless packaging is damaged. Do not use when packaging is damaged.

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