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ToxiSeal Vial Adaptors with external balloon

REF **420200** 13 mm Vial / Pressure Equalization 80 ml


REF **420210** 20 mm Vial / Pressure Equalization 80 ml

REF **420220** 28 mm Vial / Pressure Equalization 80 ml

Air Filter Pore Size 0.1 um
Carbon Filter Activated Carbon Fabric
With Valve

Packaging
20 pcs / Inner Box
80 pcs / Outer Carton

NOT MADE with NATURAL RUBBER LATEX or DEHP

780093 Rev.02, Date of issue: 2021-02-08

ToxiSeal Vial Adaptors seal the closure of the vial to which it is attached. The Air Vent assembly consists of 2 layers of filtration membranes, a hydrophobic 0.1 micron sterile filter and an activated carbon filter that absorbs drug aerosols and vapor. The external balloon collects drug aerosols and volatile organic compounds released due to un-equalized pressure in the vial during compounding.

The injection site on the ToxiSeal Vial Adaptor with external balloon and all its corresponding interface membranes exhibit a dry connection with the communicating surfaces in a fluid transfer. The use of this component and its appropriate ProSeal CSTD connecting component reduces the risk of microbial ingress for up to 168 hours or 7 days.

DIRECTIONS - Use Aseptic Technique

1. Attach ToxiSeal Vial Adaptor with external balloon firmly to Vial.
2. Prior to every access, swab top of Injection Site with 70% isopropyl alcohol (15 seconds) and allow to dry (approximately 30 secs).
Dry Time is dependent on temperature, humidity, ventilation area.
3. Attach ProSeal Injector (Syringe Adaptor) Ref 421010 to the Injection Site of ToxiSeal Vial Adaptor with external balloon. Invert Vial to withdraw fluid.
4. The external balloon allows solution to be withdrawn without injecting air. Avoid pushing syringe plunger when the vial is inverted as excessive pressure may compromise the filter membranes.
5. Discard used ToxiSeal Vial Adaptor with external balloon with vial intact in accordance to disposal procedures for biohazardous materials of your facility.

EC REP

MT Promed Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert Germany
+49-6894-581020

Epic International (Thailand) Co., Ltd.
Hemaraj Eastern Seaboard Industrial Estate
500/73 Moo 2, Tasit, Pluakdaeng, Rayong
Thailand 21140

www.epic-med.com

Indications for Use:

The ProSeal Closed System drug Transfer Device (CSTD) mechanically prohibits environmental contaminants from entering the system and the escape of drug or vapor concentrations from the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols, and spills. The ProSeal system also prevents the introduction of microbial contaminants into the drug or fluid path for up to 168 hours or 7 days.

Contraindications

The device is contraindicated whenever:

- The drug to be prepared is contraindicated to polycarbonate, ABS, stainless steel, polypropylene, PVC, silicone and polyisoprene.

Warnings

- Use accepted IV and pharmacy practice.
- The performance of the self-sealing membrane of the device is reduced after multiple perforations.
- Do not re-use to avoid contamination.

Precautions

- Do not use when caps and/or components are loose.
- The device is sterile unless packaging is damaged. Do not use when packaging is damaged.

CE 0123

STERILE EO