

Smartez™ Infusion Pumps - Soft Shell

Chemotherapy Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
484010	100 ml	2 ml/h	50 hrs
484020	270 ml	2 ml/h	135 hrs
484030	120 ml	4 ml/h	30 hrs
484040	270 ml	5 ml/h	54 hrs
484050	270 ml	10 ml/h	27 hrs
484060	50 ml	2 ml/h	25 hrs
484070	60 ml	0.5 ml/h	120 hrs
484080	240 ml	10 ml/h	24 hrs
484090	120 ml	30 ml/h	4 hrs
484100	100 ml	100 ml/h	60 min
484110	270 ml	1.5 ml/h	180 hrs
484120	60 ml	2 ml/h	30 hrs
484130	65 ml	0.5 ml/h	130 hrs
484140	240 ml	5 ml/h	48 hrs
484230	240 ml	2 ml/h	120 hrs
484240	240 ml	10 ml/h	24 hrs
484260	100 ml	5 ml/h	20 hrs

Smartez™ Infusion Pumps - Soft Shell

Chemo CS with ProSeal™ Infusion Duration

Model Number (REF)	Nominal Volume	Nominal Flow Rate	Nominal Time
484400	100 ml	2 ml/h	50 hrs
484410	270 ml	2 ml/h	135 hrs
484420	120 ml	4 ml/h	30 hrs
484430	270 ml	5 ml/h	54 hrs
484440	270 ml	10 ml/h	27 hrs
484450	50 ml	2 ml/h	25 hrs
484460	60 ml	0.5 ml/h	120 hrs
484470	240 ml	10 ml/h	24 hrs
484480	120 ml	30 ml/h	4 hrs
484490	100 ml	100 ml/h	60 hrs
484500	270 ml	1.5 ml/h	180 hrs
484510	60 ml	2 ml/h	30 hrs
484520	65 ml	0.5 ml/h	130 hrs
484530	240 ml	5 ml/h	48 hrs
484540	240 ml	2 ml/h	120 hrs
484550	240 ml	10 ml/h	24 hrs