

# SMARTeZ™ Pump

enabling what matters most



## Infusion Therapy

**SMART and Easy**



## Product Specifications and Ordering Information

**SMARTeZ™ Elastomeric Infusion Pumps** are designed to give clinicians and nurses the option of delivering pre-determined amounts of medication to the patient in a continuous and accurate manner either at the hospital or at home.

Elastomeric pump devices are intended to infuse medication/fluid for either continuous intravenous, subcutaneous or epidural infusion. Chemotherapy, antibiotics, anesthesia and pain management are the most common therapies where elastomeric pumps can be used, either in adult or pediatric patients.

For more technical information related to the use of the pump and factors affecting flow rates visit: [www.epic-med.com](http://www.epic-med.com)

### SMARTeZ™ PUMPS

REF No. ( Article No. )	Nominal Fill Volume(mL)	Nominal Flow Rate (mL/h)	Nominal Duration ( Hr )	Qty Per Carton	Color Code	
<b>Short Duration</b>						
481031	50	50	1	10	Green	
481091	100	50	2	10	Green	
481151	250	50	5	10	Green	
481041	100	100	1	10	White	
481111	200	100	2	10	White	
481121	250	100	2 hrs 30 min	10	White	
481101	400	100	4	10	White	
481131	250	125	2	10	Dark Blue	
481011	100	200	0.5	10	Turquoise	
481141	200	200	1	10	Turquoise	
481071	400	200	2	10	Turquoise	
481051	250	250	1	10	Yellow	
481081	500	250	2	10	Yellow	
481021	250	500	0.5	10	Violet	
<b>Long Duration</b>						
480151	65	0.5	130	10	Blue	
480181	100	0.5	200	10	Blue	
480191	270	1	270	10	Pink	
480121	100	1.5	67	10	Dark Blue	
480051	60	2	30	10	Yellow	
480081	100	2	50	10	Yellow	
480101	120	2	60	10	Yellow	
480161	270	2	135	10	Yellow	
480171	300	2	150	10	Yellow	
480061	120	4	30	10	Grey	
480131	270	4	68	10	Grey	
480141	400	4	100	10	Grey	
480011	60	5	12	10	Orange	
480021	80	5	16	10	Orange	
480031	125	5	25	10	Orange	
480091	270	5	54	10	Orange	
480111	400	5	80	10	Orange	
480041	270	10	27	10	Light Green	
480071	400	10	40	10	Light Green	
<b>Chemo</b>						
484011	100	2	50	10	yellow	
484021	270	2	135	10	yellow	
484031	120	4	30	10	Grey	
484041	270	5	54	10	Orange	
484051	270	10	27	10	Light Green	



[www.epic-med.com](http://www.epic-med.com)



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert Germany